


Mirante

FA/ICG/OCT



The Ultimate Multimodal Imaging Platform
State-of-the-art SLO/OCT Combo

Ultra Wide Field x Ultra HD image

A stellar combination of 163° ultra wide field x ultra 4K HD incorporated in the Mirante achieves a wider, enhanced view of the retinal structure and vasculature with unparalleled clarity. Ultra wide field image is available with the optional wide-field adapter. **Measured from the center of the eye.

FlexTrack

The new FlexTrack technology improves imaging quality.

Images courtesy of Giovanni Staurenghi, MD and Maurizio Cori, MD, Luigi Sacco Hospital, University of Milan, Italy


SLO/OCT model



www.nidek.com

Multimodal Imaging Structure & Function

Evaluate retinal structure and function simultaneously using combined OCT and Microperimetry images



OCT

Microperimeter



Mirante SLO/OCT



RS-3000 Advance 2



MP-3



www.nidek.com